



Diode laser Mantis HR Black Mamba is your best ally for laser epilation for all phototypes.

The **highest degree of performance** using a **single handpiece** that operates at two levels of depth.

This allows to treat a **wider range of clients** throughout the year, and is delicate on the skin.










All the Mantis technology in just one, single handpiece.

The patented source and low weight allow the **Black Mamba handpiece** the highest degrees of performance, making the treatment a unique, pleasurable experience for the operator.





New Dynamic Cooling System


An exclusive dynamic cooling system allows HR Black Mamba to maintain the highest degrees of performance throughout the treatment, guaranteeing the same temperature from the first pulse to the last on the skin, even in the most sensitive areas.



The thermal capacity of the laser source is enhanced by the ACS - Active Cooling System, guaranteeing virtually unlimited treatment times and increasing the thermal capacity of the entire laser system.



The software electronically controls the temperature and measures the water flow, thus making it possible to prevent the most common faults that tend to occur in less evolved laser systems.





Safe treatment on all phototypes. Fast and painless results.

The inclusion of the 1064YAG wavelength allows a lower absorption by the melanin to guarantee deeper action.

This feature permits safe and complete epilation even on darker phototypes and on tanned skins.

Mantis HR Black Mamba technology make possible to work in a safe and effective operation, using pre-set work protocols, which guide the operator in an assisted way.

Type1 Fair skin Red or light blonde hair Numerous freckles



Fair skin Blonde hair



Dark blonde or light brown hair



Brown hair



Dark brown hair




BLACK MBA

Mantis Care System

Tele-assistance is a remote technical support **system** that provides direct or indirect interaction with the equipment. Mantis technicians intervene remotely, without the need to intervene on site, guaranteeing continual and effective supervision.



Dual wavelength. Power, control, safety.





The 808 wavelength has an intense absorption of melanin and a deep action in the follicle. The 1064 YAG wavelength has a low absorption of melanin and even a deeper photoepilation, reaching beyond the bulb, also the papilla.







RIVENDITORE AUTORIZZATO



Mantis Srl

Via Artigianale Croce, 13 42035 Castelnovo ne' Monti (RE) - Italy Tel +39 0522 612080 Fax +39 0522 729796 info@mantisitalia.it









mantisitalia.it



